


The Gold Standard to Treat Atrial Flutter




AlCath Flutter Gold

AlCath Flutter LT G and AlCath Flutter Flux G eXtra are specifically designed to deliver refined control and enhanced stability. Both catheters provide lateral stabilization of the catheter tip segment to follow isthmus anatomy smoothly without jumping off or losing effective contact.


Additionally, AlCath Flutter Gold features innovative technology, such as high-density braiding and lateral stabilization of the tip segment, to provide enhanced shaft stability and offer an optimized design for treatment of atrial flutter.

Precise Performance. Reliable Results.




Reinforced Stability

Effective tip-to-tissue contact enabled by lateral stabilization of tip segment.¹



Refined Control

High-density double braiding for refined control and direct 1:1 motion transfer.




Reduced Length

Optimized for right atrial CTI ablations with a length of 95cm.


Enhanced Control

STABLE DRAGGING



Good reach and positioning even across a complex isthmus.

EFFECTIVE FORCE APPLICATION




Effective force application, even in challenging areas close to the tricuspid annulus.


Mapping a New Frontier

Non-contact 3D Imaging & Mapping with Integrated Therapy

Speak to your Acutus representative today to discover how our complete mapping and integrated therapy portfolio can support your EP laboratory.




Advanced Features. Enhanced Efficacy.



Gold-tip Ablation

Gold-tip ablation delivers superior acute efficacy—enabling shorter procedure times.²



Advanced Active Cooling*

Delivers a low-flow rate without compromising on safety or efficacy.

*Only applies to AICath Flutter Flux G eXtra, AICath Flutter LT G is a non-irrigated catheter

Specifications & Ordering Information

Model	Tip	Curve	Spacing	Shaft Diameter	Length	Order Number
AICath Flutter Flux G eXtra	3.5 mm	Black	2-5-2 mm	7 F	95 cm	380783
AICath Flutter LT G	8 mm	Black	2-5-2 mm	7 F	95 cm	368997

ACUTUS MEDICAL, INC.

2210 Faraday Ave
Suite 100
Carlsbad, CA 92008
United States
Phone: +1-442-232-6080
Acutus.com

ACUTUS MEDICAL NV

Ikaroslaan 25
1930 Zaventem
Belgium
Phone: +32-2-669-75-00
Fax: +32-2-669-75-01

CONTACT US



ACUTUS[®]
M E D I C A L

Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events, and directions for use. acutus.com/patents

Acutus Medical®, Acutus® and AcQMap® are registered trademarks of Acutus Medical, Inc. AcQBlade® is a trademark of Acutus Medical. This material is the property of Acutus Medical and its subsidiaries. Please do not copy, forward, or distribute this document. Copyright © 2021 Acutus Medical, Inc. All rights reserved.

Specifications are subject to modification, revision and improvement.

1. AICath Flutter Post-Market Observation as summarized in CER-000-009-840_A.

2. T. Lewalter, C. Weiss, C. Mewis, W. Jung, W. Haverkamp, J. Proff, W. Bauer: An optimized approach for right atrial flutter ablation: a post hoc analysis of the AURUM 8 study. J. Interv. Card Electrophysiol. (2017) 48:159-166. DOI 10.1007/s10840-016-0202-x